




5V2A USB + Type C Flexible Electric Belly Belt Massage Personalized Comfort

Our Product Introduction

For more products please visit us on electronicsmassager.com

Basic Information

- Place of Origin: Made in China
- Brand Name: GT
- Certification: CE,FCC,ROHS
- Model Number: GT-1606
- Minimum Order Quantity: 1000
- Price: \$17.8-\$18.8/unit
- Packaging Details: Carton,Customized Packaging
- Delivery Time: 5-25 work days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 10000/month



Product Specification

- Features: Adjustable, Portable, Rechargeable
- Frequency: 50-60Hz
- Color: Black
- Temperature Levels: 3
- Power Source: USB
- Weight: 0.5 Kg
- Size: One Size
- Voltage: 5V
- Highlight: belly belt massage personalized,
5v2a belly belt massage,
5v2a belly fat massage belt



More Images



Product Description

5V2A USB + type CFlexible Electric Waist Support Belt for Personalized Comfort

Product Description:

The Belly Massage Belt is designed to target your belly area, providing a deep tissue massage that can help with blood flow and reduce inflammation. With three massage modes to choose from, you can switch up the intensity to suit your needs. The three temperature levels can also help soothe sore muscles and provide a relaxing experience.

Not only does the Belly Massage Belt help with your slimming goals, but it can also be used as a heating waist belt to help with pain relief. Whether you're dealing with sore muscles after a workout or just need a little extra heat to help with menstrual cramps, the Belly Massage Belt can provide the relief you need.

Included in the package contents are the Belly Massage Belt, a USB cable for easy charging, and a user manual to help you get started. The Belly Massage Belt is easy to use and can be adjusted to fit your waist comfortably. Simply put it on, choose your massage and temperature settings, and relax while the belt does the work.

The Belly Massage Belt is a versatile tool that can be used for a variety of purposes. Whether you're looking to target your belly area for slimming purposes or need a heating waist belt for pain relief, the Belly Massage Belt has you covered. With its durable neoprene material and three massage modes and temperature levels, the Belly Massage Belt is a great addition to any home spa kit.

Features:

Product Name: Belly Massage Belt

Voltage: 5V

Material: Neoprene

Power Consumption: 5W

Model: BM-01

Color: Black

Embossed Leather Belt

Belly Dance Accessories

Back Massager Pad

Technical Parameters:

Model	BM-01
Temperature Levels	3
Power Consumption	5W
Material	Neoprene
Massage Intensity Levels	3
Weight	0.5 Kg
Frequency	50-60Hz
Power Source	USB
Package Contents	1 X Belly Massage Belt, 1 X USB Cable, 1 X User Manual
Voltage	5V

Applications:

The Belly Massage Belt is ideal for a variety of occasions and scenarios. It can be used while sitting at a desk, relaxing at home, or even during a workout. For those who suffer from back pain, this belt doubles as a Back Massager Pad, providing relief and comfort to the entire back area.

For those who are looking to improve their physique, the GT-1606 Belly Massage Belt can also be used as a Bodybuilding Waist Trainer. The unique design of the belt helps to target the core muscles, providing a more effective workout and increased results.

Additionally, the Belly Massage Belt can also be used as a Belly Dance accessory. The three temperature levels allow for a customizable experience, making it perfect for dancers who want to warm up their core muscles before a performance.

Overall, the GT-1606 Belly Massage Belt is a versatile and effective product that can be used for a variety of occasions and scenarios.

Made in China with high-quality materials, this belt is durable and long-lasting, providing relief and comfort for years to come.

Customization:

Customize your GT Belly Massage Belt with our product customization services!

Brand Name: GT

Model Number: GT-1606

Place of Origin: Made in China

Temperature Levels: 3

Model: BM-01

Features: Adjustable, Portable, Rechargeable

Material: Neoprene

Color: Black

Make your Belly Massage Belt even more personalized by adding a Back Massager Pad or wear it as Belly Dance Accessories. Keep your waist warm with the Heated Waist Belt feature.

Support and Services:

Our Belly Massage Belt product comes with comprehensive technical support and services to ensure you have the best experience possible. Our team of experts are available to answer any questions you may have about the product, its features, and how to use it effectively. Additionally, we offer a range of services including troubleshooting, repair, and maintenance to keep your Belly Massage Belt in top condition. We also provide training and resources to help you get the most out of your product.

Packing and Shipping:

Product Packaging:

The Belly Massage Belt will be packaged in a sturdy and attractive box to ensure safe delivery. The box will be labeled with the product name, features, and usage instructions. The packaging will also include a user manual and warranty information.

Shipping:

Our standard shipping method is through reliable courier services. The shipping cost will be calculated based on the destination and weight of the package. Once the order is confirmed, the product will be shipped within 24 hours. The estimated delivery time is 3-5 business days within the United States and 7-14 business days for international orders. Customers will be provided with a tracking number to track their package.

FAQ:

Q: What is the brand name of this product?

A: The brand name of this product is GT.

Q: What is the model number of this Belly Massage Belt?

A: The model number of this Belly Massage Belt is GT-1606.

Q: Where is this Belly Massage Belt made?

A: This Belly Massage Belt is made in China.

Q: Is this Belly Massage Belt suitable for all body types?

A: Yes, this Belly Massage Belt is designed to be adjustable and suitable for all body types.

Q: How often should I use this Belly Massage Belt?

A: We recommend using this Belly Massage Belt for 15-20 minutes per session, 2-3 times per week for best results.

gotole Shenzhen Zhizun Industrial Co., Ltd.

+8613828753486

szgotolesq@163.com

electronicsmassager.com

Room 301,3rd Floor, Office Building 16, Zhongpengcheng Industrial Park., No.3 Zhenxing Road, Matian Street,518107,Guangming District, Shenzhen,Gongdong. China.