### Remote Controlled Leg Massager With Heat And Adjustable Wraps

### Basic Information

· Place of Origin: Made in China

• Brand Name: GT

CE,FCC,ROHS Certification: Model Number: GT-1509 Minimum Order 1000

Quantity:

Our Product Introduction

more products please visit us on electronicsmassager.com

• Price: USD\$11.9-13.2\$/unit

 Packaging Details: Carton, Customized Packaging

• Delivery Time: 5-25 work days

• Payment Terms: T/T Western Union Money Gram

 Supply Ability: 10000/month



### **Product Specification**

· Remote Control:

16.8 X 15.3 X 7.3 Inches • Dimensions:

• Auto Shut-Off: Yes Voltage: 110 Volts • Adjustable Leg Wraps: Yes Material: ABS Model: RF-ALM070

24 Watts • Highlight: remote controlled leg massager with heat,

adjustable leg wrap massager, adjustable leg massager with heat



### More Images

• Wattage:











### **Product Description**

Remote Controlled Leg Massager with Adjustable Wraps

### **COMPANY PROFILE**

Shenzhen Guangtong precision Industry Co., LTD



# ABOUT US

Shenzhen Guangtong precision Industry Co., LTD

Shenzhen Guangtong precision Industry Co., LTD professional high-tech enterprise engaged in the operation of electronic massage products. With over 20 years of development, Guangtong now has a plant area of more than 2000 square meters. We specialize in the sales of massage products such as cupping instruments, neck and shoulder massagers, massage pillows, massage hammers, massage guns, eye beauty massagers, slimming and weight loss massage belts, massage pads, foot massagers, and more. Our products have been sold in North America, Europe, Asia, and other countries.

LEARN MORE >>

 $\bigcirc$ 

 $\subset$ 

0

 $\leq$ 

刀

T I

 $\circ$ 

0

more products please visit us on electronicsmassager.com













### **CERTIFICATE DISPLAY**



### WE ARE AN ENTERPRISE INTEGRATING R & D, PRODUCTION AND SALES.

We will continue to accelerate the development and innovation of new products, and develop and promote new markets. Amthi is willing to build a new green energy platform with friends.















I W 0 Z I 0 0 S





## F A Q

#### Q1: Could you suggest a design based on our requirements or samples?

Of course, we are professional. We can test samples free of charge according to your requirements.

#### Q2: Do you help with the operation of the product?

We have a product instruction manual. Of course, if you have any problems, you can find me through the mail or online chat window to solve the operation problems for you!

#### O3: How do I process orders?

Let us know your requirements and applications, and then we provide you with detailed quotation and logistics according to your requirements, and finally we confirm the order details together!

#### Q4: How soon can you deliver the equipment?

We will deliver your equipment within 3–5 days (there are a lot of stock), custom or mold processing needs to take 7~15 days.

#### Q5: About packaging and shipping

Provide strong plywood packing machine, safe and in line with international export standards.

#### Q6: How do you ensure we receive a good product?

Our equipment is their own procurement of raw materials, their own installation and commissioning. We understand all aspects of our equipment and will not go through a third party, and we will send you the equipment debugging video before shipment, to ensure that there will be no quality problems on your hands, causing trouble to your work!

### **Product Description:**

The Electric Scraping Massager is made of high-quality ABS material, which ensures durability and longevity. This material is known for its strength and resistance to wear and tear, making it the perfect choice for a product that will be used frequently and for extended periods.

This leg massager is designed with a sleek and modern white finish that will complement any decor. Its compact dimensions, measuring 16.8 x 15.3 x 7.3 inches, make it easy to store when not in use, while still providing ample coverage for your legs during a massage. The Electric Scraping Massager offers a range of massage etchniques to help you relax and soothe sore muscles. Two of the most popular techniques used in this product are air compression and shiatsu. With air compression, the leg wraps inflate and deflate to apply pressure to the legs, promoting circulation and reducing swelling. Shiatsu massage, on the other hand, uses rolling, kneading, and pressing motions to relieve tension and pain in the muscles.

If you're looking for a complete massage experience, you might also consider pairing the Electric Scraping Massager with a Foot Bath

If you're looking for a complete massage experience, you might also consider pairing the Electric Scraping Massager with a Foot Bath Massager. This will provide a full-body relaxation experience, as the foot bath massager will help to soothe tired feet and promote overall relaxation.

The Electric Scraping Massager can also be used in conjunction with a Body Slimmer Massager to help reduce the appearance of cellulite and tone your muscles. This combination works by breaking down fat cells and promoting circulation, resulting in smoother, firmer

Overall, the Electric Scraping Massager is a versatile and effective product that can help you enjoy the benefits of massage therapy from the comfort of your own home. With its adjustable leg wraps, durable ABS material, and range of massage techniques, this leg massager is an excellent investment for anyone looking to relieve pain, reduce swelling, and promote relaxation.

### Features:

Product Name: Leg Massager

Color: White

Material: ABS

Massage Intensity Levels: 3

Model: RF-ALM070

Auto Shut-Off: Yes

Product Type: Anti Cellulite Electric Massager

Usage: Handheld Muscle Massager

### **Technical Parameters:**

Voltage	110 Volts
Auto Shut-Off	Yes
Remote Control	Yes
Wattage	24 Watts
Color	White
Model	RF-ALM070
Adjustable Leg Wraps	Yes
Power Source	Corded Electric
Material	ABS
Massage Intensity Levels	3

### Applications:

The GT-1509 Leg Massager is perfect for a variety of occasions and scenarios. Athletes and fitness enthusiasts can use it to relieve muscle soreness and improve blood circulation after a workout. People with desk jobs that require long hours of sitting can use it to prevent leg fatigue and swelling. The massager can also be used by individuals with medical conditions such as varicose veins and

Impledema to alleviate symptoms.

The GT-1509 Leg Massager is not limited to just the legs, it can also be used on other body parts such as the arms, back, and feet. Its versatile design makes it a great addition to any home. The massager can also be used in conjunction with other massage products such as a handheld muscle massager, electric scraping massager, or foot bath massager for a complete at-home spa experience.

For businesses looking to purchase the GT-1509 Leg Massager, the minimum order quantity is 1000 units with a price range of USD \$11.9-13.2 per unit. The packaging details can be customized to fit the needs of the business. The delivery time ranges from 5-25 work days depending on the quantity ordered and the location of the business. Payment terms include T/T, Western Union, or Money Gram. The supply ability of the GT-1509 Leg Massager is 10000 units per month.

### **Customization:** Support and Services:

Our Leg Massager product comes with comprehensive technical support and services to ensure our customers get the most out of their purchase. Our team of knowledgeable product experts are available to assist with any questions or concerns regarding the operation, maintenance, or troubleshooting of the Leg Massager.

In addition, we offer a range of services to enhance the customer experience, including product demonstrations, user manuals, and

warranty support. Our goal is to provide our customers with the highest level of support and service to ensure their satisfaction with our Leg Massager product.

### Packing and Shipping:

Product Packaging:
Our Leg Massager comes in a sturdy and durable box to ensure safe delivery. The box includes:

Leg Massager unit

Power adapter

Remote control

User manual

Shipping

We offer free shipping to any location within the United States. Orders are processed within 24 hours and shipped via standard delivery which usually takes 3-5 business days. Expedited shipping is also available for an additional fee

### FAQ:

- Q: What is the brand name of this Leg Massager?
  A: The brand name of this Leg Massager is GT.
- Q: What is the model number of this Leg Massager?
- A: The model number of this Leg Massager is GT-1509.
- Q: Where is this Leg Massager made? A: This Leg Massager is made in China.
- Q: What certifications does this Leg Massager have?
- A: This Leg Massager has CE, FCC, and ROHS certifications.

  Q: What is the minimum order quantity for this Leg Massager?

- A: The minimum order quantity for this Leg Massager is 1000 units.

  Q: What is the price range for this Leg Massager?

  A: The price range for this Leg Massager is USD\$11.9-13.2\$ per unit.
- Q: What are the packaging details for this Leg Massager?
  A: The packaging details for this Leg Massager are carton and customized packaging.
- Q: What is the delivery time for this Leg Massager?
  A: The delivery time for this Leg Massager is 5-25 work days.
- Q: What are the payment terms for this Leg Massager?
- A: The payment terms for this Leg Massager are T/T, Western Union, and Money Gram. Q: What is the supply ability for this Leg Massager?
- A: The supply ability for this Leg Massager is 10000 units per month.

Gotole Shenzhen Zhizun Industrial Co., Ltd.



szgotolesq@163.com electronicsmassager.com

Room 301,3rd Floor, Office Building 16, Zhongpengcheng Industrial Park., No.3 Zhenxing Road, Matian Street,518107,Guangming District, Shenzhen,Gongdong. China.