

110Volt Compression Leg Wraps Circulation FCC Leg Massager With Heat

Basic Information

Place of Origin: Made in China

Brand Name: GT

Certification: CE,FCC,ROHSModel Number: GT-1509

 Minimum Order Quantity:

er 1000

• Price: USD\$11.9-13.2\$/unit

Packaging Details: Carton, Customized Packaging

• Delivery Time: 5-25 work days

• Payment Terms: T/T, Western Union, Money Gram

• Supply Ability: 10000/month



Care for health starts with legs

Product Specification

• Massage Techniques: Air Compression, Shiatsu

Color: WhiteAdjustable Leg Wraps: YesAuto Shut-Off: Yes

Power Source: Corded ElectricModel: RF-ALM070

Remote Control: YesVoltage: 110 Volts

Highlight: 110volt compression leg wraps circulation,

fcc leg massager with heat,

compression leg wraps circulation fcc



More Images











Elegant life Healt

Product Description

110 Volt Air Compression and Shiatsu Air Compression Leg Massager with Adjustable Leg Wraps

Product Description:

With its handheld body massager design, you can easily control the intensity and mode of your massage using the remote control. The RF-ALM070 model offers three different massage intensity levels and modes, giving you the flexibility to choose the perfect massage for your specific needs

Whether you are looking for a gentle massage to help you relax or a deeper tissue massage to help relieve soreness and tension, this massager has you covered. The adjustable leg wraps ensure that you can customize your massage experience to suit your individual

The RF-ALM070 Anti-Cellulite Electric Leg Massager is perfect for anyone who spends long hours on their feet or who suffers from poor circulation. The three different massage modes allow you to target specific areas of your legs, helping to improve blood flow and reduce the appearance of cellulite.

With its sleek and modern design, this massager is easy to use and store when not in use. So why wait? Treat yourself to the ultimate leg massage experience with the RF-ALM070 Anti-Cellulite Electric Leg Massager today!

Features:

Product Name: Leg Massager

Wattage: 24 Watts Color: White

Adjustable Leg Wraps: Yes

Massage Techniques: Air Compression, Shiatsu

Dimensions: 16.8 X 15.3 X 7.3 Inches

This Electric Scraping Massager offers air compression and shiatsu massage techniques to help relieve tension and relax your legs. With adjustable leg wraps, the Leg Massager fits comfortably on any leg size. It also comes in a sleek white color and has a wattage of 24 Watts. Additionally, this Anti Cellulite Electric Massager provides heat to your feet for added comfort.

Technical Parameters:

Model	RF-ALM070
Wattage	24 Watts
Massage Techniques	Air Compression, Shiatsu
Auto Shut-Off	Yes
Adjustable Leg Wraps	Yes
Color	White
Material	ABS
Remote Control	Yes
Voltage	110 Volts
Dimensions	16.8 X 15.3 X 7.3 Inches

Applications:

The GT Leg Massager GT-1509 is a versatile product that can be used in a variety of occasions and scenarios. Below are some product

After a Long Day at Work

After a long and tiring day at work, your legs can become sore and swollen. The GT Leg Massager GT-1509 is the perfect solution to relieve your tired legs. Simply wrap the adjustable leg wraps around your legs and select your desired massage mode using the remote control. The massager will provide a soothing massage that will help to reduce swelling and relieve muscle tension. **During Pregnancy**

 $Pregnancy \ can be \ a \ challenging \ time \ for \ women, \ and \ leg \ cramps \ and \ swelling \ are \ common \ symptoms. \ The \ GT \ Leg \ Massager \ GT-1509$ can be used during pregnancy to provide relief from these symptoms. The adjustable leg wraps can be adjusted to fit your growing belly, and the massager will provide a gentle massage that will help to reduce swelling and relieve muscle tension. After Exercise

After a strenuous workout, your legs can become sore and stiff. The GT Leg Massager GT-1509 is the perfect product to help you recover after exercise. The massager can be used to provide a deep tissue massage that will help to relieve muscle soreness and speed up recovery time.

Foot Bath Massager

The GT Leg Massager GT-1509 can also be used as a foot bath massager. Simply fill a bowl with warm water and place the leg wraps in the water. Select your desired massage mode using the remote control, and the massager will provide a soothing massage that will help to relieve foot pain and improve circulation.

The GT Leg Massager GT-1509 has a wattage of 24 watts and offers 3 massage modes. The adjustable leg wraps ensure that the massager can be used by people of all sizes. The product is available in white color and comes with a remote control for easy operation. The GT Leg Massager GT-1509 is supplied by GT, and the company has a supply ability of 10000/month. The payment terms include T/T, Western Union, and Money Gram.

Customization:

- Brand Name: GT
- Model Number: GT-1509
- Place of Origin: Made in China
- Certification: CE, FCC, ROHS
- Minimum Order Quantity: 1000
- Price: USD\$11.9-13.2\$/unit
- Packaging Details: Carton, Customized Packaging
- Delivery Time: 5-25 work days
- Payment Terms: T/T, Western Union, Money Gram
- Supply Ability: 10000/month Auto Shut-Off: Yes
- Power Source: Corded Electric
- Adjustable Leg Wraps: Yes
- Color: White
- Massage Modes: 3

Make your Handheld Muscle Massager, Electric Massager, or Foot Bath Massager even better with our Product Customization Services!

Our Leg Massager product comes with comprehensive technical support and services to ensure you have the best experience possible. Our team of experts is available to answer any questions you may have about the product, including its features, functions, and operation. We also offer troubleshooting assistance, should you encounter any issues with your Leg Massager. Our team can help you identify and resolve any problems you may be experiencing, ensuring that your product is working optimally at all times.

In addition to technical support, we also offer a range of services to help you get the most out of your Leg Massager. These include product training, maintenance and repair services, and replacement parts and accessories

With our comprehensive support and services, you can feel confident in your purchase of the Leg Massager, knowing that you have access to the resources you need to keep it running smoothly and effectively

Packing and Shipping:

Product Packaging:

The Leg Massager will be securely packed in a sturdy cardboard box to ensure safe delivery. The product will be wrapped in bubble wrap to prevent any damage during shipping.

Shipping:

We offer free standard shipping for all orders within the United States. The Leg Massager will be shipped within 2-3 business days of order confirmation and will arrive within 5-7 business days. Expedited shipping options are available at an additional cost.

FAQ:

- Q: What is the brand name of this leg massager?
- A: The brand name of this leg massager is GT.
- Q: What is the model number of this leg massager?
- A: The model number of this leg massager is GT-1509.
- Q: Where is this leg massager made?
- A: This leg massager is made in China.
- Q: What certifications does this leg massager have?
 A: This leg massager has CE, FCC, and ROHS certifications.
- Q: What is the minimum order quantity for this leg massager?
- A: The minimum order quantity for this leg massager is 1000 units.
- Q: What is the price range for this leg massager?
- A: The price range for this leg massager is USD\$11.9-13.2\$/unit.
- Q: What are the packaging details for this leg massager?
 A: The packaging details for this leg massager are carton and customized packaging.
- Q: What is the delivery time for this leg massager?
- A: The delivery time for this leg massager is 5-25 work days.
- Q: What are the payment terms for this leg massager?

 A: The payment terms for this leg massager are T/T, Western Union, and Money Gram.
- Q: What is the supply ability for this leg massager?
- A: The supply ability for this leg massager is 10000 units per month.

Gotole Shenzhen Zhizun Industrial Co., Ltd.



+8613828753486



szgotolesq@163.com



electronicsmassager.com

Room 301,3rd Floor, Office Building 16, Zhongpengcheng Industrial Park., No.3 Zhenxing Road, Matian Street,518107, Guangming District, Shenzhen, Gongdong. China.